



**COLLEGE OF SCIENCE, TECHNOLOGY & APPLIED ARTS**

**OF TRINIDAD & TOBAGO**

**School of Nursing Health and Environmental Sciences  
Department of Natural & Life Sciences**


---

## **CHEM 092 – Introduction to Concepts in Chemistry II**

---

**Student Manual**

Edited: 29/02/2016



## A Word of Welcome

The School of Nursing Health and Environmental Sciences and the Department of Natural & Life Sciences welcome you to the practical component of the course. The experiments that form part of this course have been selected to:

- Reinforce your knowledge and understanding of the theoretical content of your first year Chemistry course.
- Develop skills such as observing, recording, measuring and reporting of results.
- Teach you special techniques, such as weighing, pipetting, titrating, filtration etc, so that you will become familiar with the basics of laboratory work.



## TABLE OF CONTENTS

A WORD OF WELCOME .....	2
TABLE OF CONTENTS .....	3
INTRODUCTION .....	4
SAFETY IN THE CHEMISTRY LABORATORY.....	7
WASTE DISPOSAL.....	13
MANAGING YOUR TIME DURING THE PRACTICAL.....	14
GETTING STARTED.....	15
LAB 1: TITRATION METHOD.....	18
LAB 2: DETERMINATION OF CONCENTRATION OF AN UNKNOWN SOLUTION BY TITRATION METHOD.....	20
LAB 3: REDOX TITRATION .....	22

## Introduction

### General

Always read and study the relevant experiment before coming to the laboratory. Check out appropriate references to help you understand the procedures you will be carrying out.

### Attendance at all laboratory sessions is compulsory.

Learning lab techniques must come from *experiencing* lab techniques. Experience is gained from carrying out experiments. You will need to plan ahead for your lab session, and think about what you are doing (and why) while you are in lab.

### Use past experiences to develop new skills.

You will gain new skills with new equipment and learn more about the chemistry theories and concepts you hear about in lectures, but the skills that we expect you to develop the most during these experiments are general scientific skills that should be transferable to any other scientific discipline and to your daily life. These skills developed in the laboratory are fundamental to the world of work because industries depend on the laboratory to maintain a level of quality, which enhances the company's reputation, as well as for product development.

The lab is a part of the learning process.

### LABORATORY SCHEDULE:

The Lab schedules are posted at the start of the semester on notice boards, lab tech stations and on the lab doors. Please refer to them to find out the dates and times of your lab sessions.

Once the semester begins, you will have to balance your laboratory study time with all of the other demands that you have. Most of the experiments have been scheduled so that the average student can finish everything (including the report) required in the laboratory period.

**You are expected to attend all scheduled laboratory periods and make-up labs will not be scheduled.**

# SONHES Multi-Science Laboratory Rules

---

## Before Entering the Lab

**Proper attire is required for entry into the lab\***

- Lab Coat
  - Long sleeve (wrist length)
  - Knee Length
  - Drill cotton
- Long pants or long skirt (ankle length)
- Shoes that fully cover your feet
- Long hair tied back
- No dangling jewelry or accessories

## Punctuality

- Students are not allowed entry into lab 15 minutes **after** the start of the lab session.


## During the Lab Session

- No eating, No drinking, No chewing of gum, No smoking
- Mobile phones and electronic devices (laptops, tablets etc.) are not to be used in the lab.
- Students are expected to follow the lab rules and to follow the guidance of the lab technicians while in the lab.
- Students are expected to have their own lab manuals at each lab session. No new manuals will be issued to students who forgot to bring theirs.

## Leaving the Lab

- Students are required to wash all glassware, return items to student cupboards as well as dispose of gloves, tissue, scrap paper etc and wipe countertops before leaving the lab.
- Lab coats are to be taken off when leaving the lab (including trips to the bathroom).

**\*LAB LECTURERS OR LAB TECHNICIANS MAY PREVENT ENTRY OR ASK STUDENTS TO LEAVE THE LAB IN THE EVENT OF STUDENT NON-COMPLIANCE TO THE LABORATORY RULES.**



- ✓ Wrist length sleeves for coverage
- ✓ Protects clothes
- ✓ Able to button at front
- ✓ Thick material for protection from spills
- ✓ Durable, long lasting protection for future labs
- ✓ Feel like a scientist!!



- ✗ Minimal barrier against chemicals and microbes
- ✗ Skin exposed
- ✗ Thin material allows chemicals to permeate and reach skin
- ✗ Hoodies can transmit contaminants to friends & family




- ✓ Feet fully covered!
- ✓ Comfortable
- ✓ Ankle length Pants
- ✓ Protection from spills and broken glassware



- ✗ Feet exposed
- ✗ Limited protection from spills
- ✗ Possibility of cuts from broken glassware
- ✗ Possible exposure to pathogens

## SAFETY IN THE CHEMISTRY LABORATORY

The laboratory is a hazardous environment in which to work. The hazards are often unavoidable since scientists regularly have to use hazardous materials and equipment. However, with sensible precautions the laboratory is probably no more dangerous than your home. You **MUST AT ALL TIMES OBSERVE ALL RULES** stated in this manual.



## Proper Lab wear

Gloves, lab coat, proper lab shoes

Barrier to prevent contact of contaminants with skin

## Be Punctual

Check the schedule regularly, arrive to lab on time

Be present for all safety information regarding the lab

## Good Hygiene

Wash hands when leaving the lab, swab counter tops after lab

Prevents cross contamination to surfaces outside of lab

## Proper use of Equipment

Check electric cords before use

Always use equipment for its intended purpose

## Essential Rules for Laboratory Safety

The essential rules for laboratory safety can be expressed under two simple headings:

### **ALWAYS AND NEVER**

#### **ALWAYS**

- ♦ Familiarize yourself with the laboratory safety procedures.
- ♦ Dress appropriately.
- ♦ Wash your hands before leaving the laboratory.
- ♦ Read the instructions for the laboratory protocol **carefully** before starting any experiment.
- ♦ Handle all chemicals with great care.
- ♦ Keep your working area tidy.
- ♦ Immediately inform your instructor about spills.

#### **NEVER**

- ♦ Eat or drink in the laboratory. CHEWING OF GUM IS STRICTLY PROHIBITED!
- ♦ Smoke in the laboratory
- ♦ Inhale, taste or sniff chemicals.
- ♦ Fool around or distract neighbours.
- ♦ Carry out unauthorized experiments.

## Laboratory Safety Procedures

Your laboratory will have certain procedures, which you must be familiar with. Make sure you know where all the exits from the laboratory are, in the event of an evacuation because of fire or other incident. Locate the SHOWERS and eyewash stations and know the type the fire extinguishers that are in the laboratory and how to operate them.

## Dress Code

- Wear a Lab coat.
- Wear shoes that cover your entire feet (sneakers) Avoid open footwear & high heels. So sandals, slipper and ballet flats are unacceptable.
- Wear approved safety goggles or safety glasses when in the lab.
- Avoid floppy garments; avoid things that dangle. These get tangled up in equipment or glassware and cause accidents.
- Avoid long, loose hair styles for the same reason. When Bunsen burners are in use, long hair sometimes catches on fire.

## Laboratory Techniques

### Handling of Acids or Bases

- In addition to the dress code above, personal protective equipment, such as Vinyl" gloves or "nitrile" GLOVES SHOULD be worn.
- For even more dangerous liquids, thick "nitrile" gloves should be used. These are very safe.
- **DO NOT add WATER to CONCENTRATED ACID.** The heat generated will cause splattering. If necessary to prepare certain solutions, **DO** add ACID to WATER (instead of the reverse order of addition). The heat generated will be less, but splattering still may occur. A good practice in all lab operations is to keep things at arm's length.

### Become acquainted with handling laboratory safety equipment e.g.

- Fire Extinguishers
- Fire Blanket
- Eye-wash Fountain
- First-Aid Kit

### Become acquainted with using of the FUME HOOD

- Use the HOOD for reactions that give off vapours, especially smelly vapours.
- The draft of the HOOD will sweep away vapours so that the lab itself maintains reasonable air quality.

### Locate the Safety Shower and understand how it works

- Shower should be used for ***dire EMERGENCY*** only!
- If you are ***ON FIRE***, or suffer a ***massive spill*** of a **dangerous chemical**, and need to get it off rapidly.

**Position yourself under the safety shower and pull the handle—a deluge of water will result.**

### Immediately clean up Broken Glass

- Sweep it up right away
- Place the broken glass in a "SHARP'S CONTAINER."
- This is a thick walled carton that will be sealed and discarded as such.


### Student Use of the Laboratory

- No UNAUTHORIZED experiments! These are terrible dangers in unskilled hands
- No EATING or DRINKING in the LAB. A good practice is to assume everything in the lab is toxic.
- DO NOT drink from lab equipment. NOT EVEN distilled water.
- When using pipettes, **DO NOT** suck the liquid into the pipette directly by mouth.
- Frequent mistakes, even by skilled lab workers, lead to the liquid reaching the mouth.
- Do NOT smell the flask directly.

### Housekeeping:


- **STUDENTS MUST CLEAN UP AFTER THEMSELVES.**
- Paper and trash must not be left in the room, thrown into drawers, sinks, on the floor, etc.
- When the equipment boxes are on carts, students who take a box from the cart to do an experiment must put everything back into the box neatly and put the box back on the cart.

If you absolutely must test the odour, carefully waft the vapour from the flask toward your nose with your hand, keeping the flask quite distant from your face.



***Correct way of wafting odors.***

Make sure that reaction tubes, e.g. test tubes, are not directed toward yourself or other persons. The chemicals may splatter out the tube.




## WASTE DISPOSAL

We all have a responsibility towards a clean and safe environment. The Chemistry Department's policy dealing with the disposal of laboratory waste is in accordance with the Occupational Health and Safety Act, which is a law of our country. We do not expect you to understand every aspect of this law, but as responsible citizen and laboratory user you must be aware of the correct procedures for disposal of the different classes of laboratory waste. You will be given specific instructions on how to dispose of your chemical wastes during each pre-lab briefing. Here are also some general rules that you must take note of:

1. Place broken glass into the specially marked bins.
2. Drain all harmless chemicals in solution into the sink followed by plenty of running water.
3. Pour toxic inorganic waste solutions into the specially labelled container in the fume cupboard.
4. Discard paper and any other solid waste into the bin.
5. Ensure that matches are extinguished before disposing of them in this way.
6. Shut off all gas and water lines when not in use.

## MANAGING YOUR TIME DURING THE PRACTICAL

1. It is recommended that you arrange your workbench according to the diagram below. Keeping your workspace organized and free of clutter will save you time and frustration, and contributes greatly to laboratory safety.



2. When collecting chemicals:
  - Choose the size of your container according to the volume of chemical that you will be collecting, for instance: for 100 cm<sup>3</sup> of solution, use a 250 cm<sup>3</sup> beaker and for 10 cm<sup>3</sup> use a test tube.
  - Before collecting chemicals and/or solutions from the dispensary, mark each container with the name or formula of the chemical to be collected. This will prevent any mix-ups later on.
3. Use a small notebook for jotting down masses, measurements and observations. Writing these bits of information on slips of scrap paper that can easily be lost is not only unprofessional, but also risky because it means that your whole afternoon's practical work can be wasted because of lost data.

## GETTING STARTED

For each laboratory period, including the first, you will be required to:

- ◆ Study the entire experimental outline in the **Laboratory Manual**. It includes specific directions concerning laboratory philosophy and protocols.
- ◆ Prepare your **Laboratory Notebook**. Include a summary and a procedural outline for your experiment (see "Rules for Keeping Your Laboratory Notebook" and Instruction in the experiment outline).
- ◆ Bring your safety goggles or use the safety goggles provided. These offer adequate protection against accidental splashing of corrosive chemicals.
- ◆ Come to lab dressed appropriately

### A GOOD PRACTICE:


Read the experimental procedure ahead of lab.

Avoid horseplay. In a laboratory setting, horseplay, even if good-natured, is absolutely unacceptable.

- **No pushing!**
- **No shoving!**

At the end of the lab period: Exit the lab in an orderly manner. Again: no running, no pushing, AND NO shoving.

Figure 1: Apparatus that you will be using in the laboratory



Atomic number																																																																																																																																																																																																																																																																																																																																																																																																																			
Symbol																																																																																																																																																																																																																																																																																																																																																																																																																			
Atomic weight																																																																																																																																																																																																																																																																																																																																																																																																																			
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18																																																																																																																																																																																																																																																																																																																																																																																																		
1 <b>H</b> 1.008	4 <b>Be</b> 9.012	9 <b>Li</b> 6.941	12 <b>Mg</b> 24.31	19 <b>K</b> 39.10	20 <b>Ca</b> 40.08	21 <b>Sc</b> 44.96	22 <b>Ti</b> 47.88	23 <b>V</b> 50.94	24 <b>Cr</b> 52.00	25 <b>Mn</b> 54.94	26 <b>Fe</b> 55.85	27 <b>Co</b> 58.93	28 <b>Ni</b> 58.69	29 <b>Cu</b> 63.55	30 <b>Zn</b> 65.39	31 <b>Ga</b> 69.72	32 <b>Ge</b> 72.61	33 <b>As</b> 74.92	34 <b>Se</b> 78.96	35 <b>Br</b> 79.90	36 <b>Kr</b> 83.80	37 <b>Rb</b> 85.47	38 <b>Sr</b> 87.62	39 <b>Y</b> 88.91	40 <b>Zr</b> 91.22	41 <b>Nb</b> 92.91	42 <b>Mo</b> 95.94	43 <b>Tc</b> 98.91	44 <b>Ru</b> 101.1	45 <b>Rh</b> 102.9	46 <b>Pd</b> 106.4	47 <b>Ag</b> 107.9	48 <b>Cd</b> 112.4	49 <b>In</b> 114.8	50 <b>Sn</b> 118.7	51 <b>Sb</b> 121.8	52 <b>Te</b> 127.6	53 <b>I</b> 126.9	54 <b>Xe</b> 131.3	55 <b>Cs</b> 132.9	56 <b>Ba</b> 137.3	57 <b>Fr</b> 223.0	58 <b>Ra</b> 226.0																																																																																																																																																																																																																																																																																																																																																																								
6 <b>C</b> 12.01	11 <b>B</b> 10.81	14 <b>Al</b> 26.98	15 <b>Si</b> 28.09	16 <b>P</b> 30.97	32 <b>Ge</b> 72.61	33 <b>As</b> 74.92	34 <b>Se</b> 78.96	35 <b>Br</b> 79.90	36 <b>Kr</b> 83.80	37 <b>Rb</b> 85.47	38 <b>Sr</b> 87.62	39 <b>Y</b> 88.91	40 <b>Zr</b> 91.22	41 <b>Nb</b> 92.91	42 <b>Mo</b> 95.94	43 <b>Tc</b> 98.91	44 <b>Ru</b> 101.1	45 <b>Rh</b> 102.9	46 <b>Pd</b> 106.4	47 <b>Ag</b> 107.9	48 <b>Cd</b> 112.4	49 <b>In</b> 114.8	50 <b>Sn</b> 118.7	51 <b>Sb</b> 121.8	52 <b>Te</b> 127.6	53 <b>I</b> 126.9	54 <b>Xe</b> 131.3	55 <b>Cs</b> 132.9	56 <b>Ba</b> 137.3	57 <b>Fr</b> 223.0	58 <b>Ra</b> 226.0																																																																																																																																																																																																																																																																																																																																																																																				
7 <b>N</b> 14.01	8 <b>O</b> 16.00	17 <b>Cl</b> 35.45	18 <b>Ar</b> 39.95	35 <b>Br</b> 79.90	36 <b>Kr</b> 83.80	37 <b>Rb</b> 85.47	38 <b>Sr</b> 87.62	39 <b>Y</b> 88.91	40 <b>Zr</b> 91.22	41 <b>Nb</b> 92.91	42 <b>Mo</b> 95.94	43 <b>Tc</b> 98.91	44 <b>Ru</b> 101.1	45 <b>Rh</b> 102.9	46 <b>Pd</b> 106.4	47 <b>Ag</b> 107.9	48 <b>Cd</b> 112.4	49 <b>In</b> 114.8	50 <b>Sn</b> 118.7	51 <b>Sb</b> 121.8	52 <b>Te</b> 127.6	53 <b>I</b> 126.9	54 <b>Xe</b> 131.3	55 <b>Cs</b> 132.9	56 <b>Ba</b> 137.3	57 <b>Fr</b> 223.0	58 <b>Ra</b> 226.0	59 <b>La</b> 138.9	60 <b>Ce</b> 140.1	61 <b>Pr</b> 140.9	62 <b>Sm</b> 150.4	63 <b>Eu</b> 152.0	64 <b>Gd</b> 157.3	65 <b>Tb</b> 158.9	66 <b>Dy</b> 162.5	67 <b>Ho</b> 164.9	68 <b>Er</b> 167.3	69 <b>Tm</b> 168.9	70 <b>Yb</b> 173.0	71 <b>Ac</b> 227.0	72 <b>Th</b> 232.0	73 <b>Pa</b> 231.0	74 <b>U</b> 238.0	75 <b>Np</b> 237.0	76 <b>Pu</b> 244.1	77 <b>Am</b> 243.1	78 <b>Cm</b> 247.1	79 <b>Bk</b> 247.1	80 <b>Cf</b> 251.1	81 <b>Es</b> 252.0	82 <b>Fm</b> 257.1	83 <b>Md</b> 258.1	84 <b>No</b> 259.1	85 <b>Uuo</b> 293	86 <b>Rn</b> 222.0	87 <b>At</b> 210.0	88 <b>Po</b> 209.0	89 <b>Bi</b> 209.0	90 <b>Pb</b> 207.2	91 <b>Tl</b> 204.4	92 <b>Hg</b> 200.6	93 <b>Uut</b> 289	94 <b>Uuq</b> 289	95 <b>Uup</b> 289	96 <b>Uuh</b> 289	97 <b>Uus</b> 289	98 <b>Uub</b> 277	99 <b>Uub</b> 272	100 <b>Uun</b> 269	101 <b>Uun</b> 269	102 <b>Uun</b> 268	103 <b>Uun</b> 268	104 <b>Uun</b> 268	105 <b>Uun</b> 268	106 <b>Uun</b> 268	107 <b>Uun</b> 268	108 <b>Uun</b> 268	109 <b>Uun</b> 268	110 <b>Uun</b> 268	111 <b>Uun</b> 268	112 <b>Uun</b> 268	113 <b>Uun</b> 268	114 <b>Uun</b> 268	115 <b>Uun</b> 268	116 <b>Uun</b> 268	117 <b>Uun</b> 268	118 <b>Uun</b> 268	119 <b>Uun</b> 268	120 <b>Uun</b> 268	121 <b>Uun</b> 268	122 <b>Uun</b> 268	123 <b>Uun</b> 268	124 <b>Uun</b> 268	125 <b>Uun</b> 268	126 <b>Uun</b> 268	127 <b>Uun</b> 268	128 <b>Uun</b> 268	129 <b>Uun</b> 268	130 <b>Uun</b> 268	131 <b>Uun</b> 268	132 <b>Uun</b> 268	133 <b>Uun</b> 268	134 <b>Uun</b> 268	135 <b>Uun</b> 268	136 <b>Uun</b> 268	137 <b>Uun</b> 268	138 <b>Uun</b> 268	139 <b>Uun</b> 268	140 <b>Uun</b> 268	141 <b>Uun</b> 268	142 <b>Uun</b> 268	143 <b>Uun</b> 268	144 <b>Uun</b> 268	145 <b>Uun</b> 268	146 <b>Uun</b> 268	147 <b>Uun</b> 268	148 <b>Uun</b> 268	149 <b>Uun</b> 268	150 <b>Uun</b> 268	151 <b>Uun</b> 268	152 <b>Uun</b> 268	153 <b>Uun</b> 268	154 <b>Uun</b> 268	155 <b>Uun</b> 268	156 <b>Uun</b> 268	157 <b>Uun</b> 268	158 <b>Uun</b> 268	159 <b>Uun</b> 268	160 <b>Uun</b> 268	161 <b>Uun</b> 268	162 <b>Uun</b> 268	163 <b>Uun</b> 268	164 <b>Uun</b> 268	165 <b>Uun</b> 268	166 <b>Uun</b> 268	167 <b>Uun</b> 268	168 <b>Uun</b> 268	169 <b>Uun</b> 268	170 <b>Uun</b> 268	171 <b>Uun</b> 268	172 <b>Uun</b> 268	173 <b>Uun</b> 268	174 <b>Uun</b> 268	175 <b>Uun</b> 268	176 <b>Uun</b> 268	177 <b>Uun</b> 268	178 <b>Uun</b> 268	179 <b>Uun</b> 268	180 <b>Uun</b> 268	181 <b>Uun</b> 268	182 <b>Uun</b> 268	183 <b>Uun</b> 268	184 <b>Uun</b> 268	185 <b>Uun</b> 268	186 <b>Uun</b> 268	187 <b>Uun</b> 268	188 <b>Uun</b> 268	189 <b>Uun</b> 268	190 <b>Uun</b> 268	191 <b>Uun</b> 268	192 <b>Uun</b> 268	193 <b>Uun</b> 268	194 <b>Uun</b> 268	195 <b>Uun</b> 268	196 <b>Uun</b> 268	197 <b>Uun</b> 268	198 <b>Uun</b> 268	199 <b>Uun</b> 268	200 <b>Uun</b> 268	201 <b>Uun</b> 268	202 <b>Uun</b> 268	203 <b>Uun</b> 268	204 <b>Uun</b> 268	205 <b>Uun</b> 268	206 <b>Uun</b> 268	207 <b>Uun</b> 268	208 <b>Uun</b> 268	209 <b>Uun</b> 268	210 <b>Uun</b> 268	211 <b>Uun</b> 268	212 <b>Uun</b> 268	213 <b>Uun</b> 268	214 <b>Uun</b> 268	215 <b>Uun</b> 268	216 <b>Uun</b> 268	217 <b>Uun</b> 268	218 <b>Uun</b> 268	219 <b>Uun</b> 268	220 <b>Uun</b> 268	221 <b>Uun</b> 268	222 <b>Uun</b> 268	223 <b>Uun</b> 268	224 <b>Uun</b> 268	225 <b>Uun</b> 268	226 <b>Uun</b> 268	227 <b>Uun</b> 268	228 <b>Uun</b> 268	229 <b>Uun</b> 268	230 <b>Uun</b> 268	231 <b>Uun</b> 268	232 <b>Uun</b> 268	233 <b>Uun</b> 268	234 <b>Uun</b> 268	235 <b>Uun</b> 268	236 <b>Uun</b> 268	237 <b>Uun</b> 268	238 <b>Uun</b> 268	239 <b>Uun</b> 268	240 <b>Uun</b> 268	241 <b>Uun</b> 268	242 <b>Uun</b> 268	243 <b>Uun</b> 268	244 <b>Uun</b> 268	245 <b>Uun</b> 268	246 <b>Uun</b> 268	247 <b>Uun</b> 268	248 <b>Uun</b> 268	249 <b>Uun</b> 268	250 <b>Uun</b> 268	251 <b>Uun</b> 268	252 <b>Uun</b> 268	253 <b>Uun</b> 268	254 <b>Uun</b> 268	255 <b>Uun</b> 268	256 <b>Uun</b> 268	257 <b>Uun</b> 268	258 <b>Uun</b> 268	259 <b>Uun</b> 268	260 <b>Uun</b> 268	261 <b>Uun</b> 268	262 <b>Uun</b> 268	263 <b>Uun</b> 268	264 <b>Uun</b> 268	265 <b>Uun</b> 268	266 <b>Uun</b> 268	267 <b>Uun</b> 268	268 <b>Uun</b> 268	269 <b>Uun</b> 268	270 <b>Uun</b> 268	271 <b>Uun</b> 268	272 <b>Uun</b> 268	273 <b>Uun</b> 268	274 <b>Uun</b> 268	275 <b>Uun</b> 268	276 <b>Uun</b> 268	277 <b>Uun</b> 268	278 <b>Uun</b> 268	279 <b>Uun</b> 268	280 <b>Uun</b> 268	281 <b>Uun</b> 268	282 <b>Uun</b> 268	283 <b>Uun</b> 268	284 <b>Uun</b> 268	285 <b>Uun</b> 268	286 <b>Uun</b> 268	287 <b>Uun</b> 268	288 <b>Uun</b> 268	289 <b>Uun</b> 268	290 <b>Uun</b> 268	291 <b>Uun</b> 268	292 <b>Uun</b> 268	293 <b>Uun</b> 268	294 <b>Uun</b> 268	295 <b>Uun</b> 268	296 <b>Uun</b> 268	297 <b>Uun</b> 268	298 <b>Uun</b> 268	299 <b>Uun</b> 268	300 <b>Uun</b> 268	301 <b>Uun</b> 268	302 <b>Uun</b> 268	303 <b>Uun</b> 268	304 <b>Uun</b> 268	305 <b>Uun</b> 268	306 <b>Uun</b> 268	307 <b>Uun</b> 268	308 <b>Uun</b> 268	309 <b>Uun</b> 268	310 <b>Uun</b> 268	311 <b>Uun</b> 268	312 <b>Uun</b> 268	313 <b>Uun</b> 268	314 <b>Uun</b> 268	315 <b>Uun</b> 268	316 <b>Uun</b> 268	317 <b>Uun</b> 268	318 <b>Uun</b> 268	319 <b>Uun</b> 268	320 <b>Uun</b> 268	321 <b>Uun</b> 268	322 <b>Uun</b> 268	323 <b>Uun</b> 268	324 <b>Uun</b> 268	325 <b>Uun</b> 268	326 <b>Uun</b> 268	327 <b>Uun</b> 268	328 <b>Uun</b> 268	329 <b>Uun</b> 268	330 <b>Uun</b> 268	331 <b>Uun</b> 268	332 <b>Uun</b> 268	333 <b>Uun</b> 268	334 <b>Uun</b> 268	335 <b>Uun</b> 268	336 <b>Uun</b> 268	337 <b>Uun</b> 268	338 <b>Uun</b> 268	339 <b>Uun</b> 268	340 <b>Uun</b> 268	341 <b>Uun</b> 268	342 <b>Uun</b> 268	343 <b>Uun</b> 268	344 <b>Uun</b> 268	345 <b>Uun</b> 268	346 <b>Uun</b> 268	347 <b>Uun</b> 268	348 <b>Uun</b> 268	349 <b>Uun</b> 268	350 <b>Uun</b> 268	351 <b>Uun</b> 268	352 <b>Uun</b> 268	353 <b>Uun</b> 268	354 <b>Uun</b> 268	355 <b>Uun</b> 268	356 <b>Uun</b> 268	357 <b>Uun</b> 268	358 <b>Uun</b> 268	359 <b>Uun</b> 268	360 <b>Uun</b> 268	361 <b>Uun</b> 268	362 <b>Uun</b> 268	363 <b>Uun</b> 268	364 <b>Uun</b> 268	365 <b>Uun</b> 268	366 <b>Uun</b> 268	367 <b>Uun</b> 268	368 <b>Uun</b> 268	369 <b>Uun</b> 268	370 <b>Uun</b> 268	371 <b>Uun</b> 268	372 <b>Uun</b> 268	373 <b>Uun</b> 268	374 <b>Uun</b> 268	375 <b>Uun</b> 268	376 <b>Uun</b> 268	377 <b>Uun</b> 268	378 <b>Uun</b> 268	379 <b>Uun</b> 268	380 <b>Uun</b> 268	381 <b>Uun</b> 268	382 <b>Uun</b> 268	383 <b>Uun</b> 268	384 <b>Uun</b> 268	385 <b>Uun</b> 268	386 <b>Uun</b> 268	387 <b>Uun</b> 268	388 <b>Uun</b> 268	389 <b>Uun</b> 268	390 <b>Uun</b> 268	391 <b>Uun</b> 268	392 <b>Uun</b> 268	393 <b>Uun</b> 268	394 <b>Uun</b> 268	395 <b>Uun</b> 268	396 <b>Uun</b> 268	397 <b>Uun</b> 268	398 <b>Uun</b> 268	399 <b>Uun</b> 268	400 <b>Uun</b> 268	401 <b>Uun</b> 268	402 <b>Uun</b> 268	403 <b>Uun</b> 268	404 <b>Uun</b> 268	405 <b>Uun</b> 268	406 <b>Uun</b> 268	407 <b>Uun</b> 268	408 <b>Uun</b> 268	409 <b>Uun</b> 268	410 <b>Uun</b> 268	411 <b>Uun</b> 268	412 <b>Uun</b> 268	413 <b>Uun</b> 268	414 <b>Uun</b> 268	415 <b>Uun</b> 268	416 <b>Uun</b> 268	417 <b>Uun</b> 268	418 <b>Uun</b> 268	419 <b>Uun</b> 268	420 <b>Uun</b> 268	421 <b>Uun</b> 268	422 <b>Uun</b> 268	423 <b>Uun</b> 268	424 <b>Uun</b> 268	425 <b>Uun</b> 268	426 <b>Uun</b> 268	427 <b>Uun</b> 268	428 <b>Uun</b> 268	429 <b>Uun</b> 268	430 <b>Uun</b> 268	431 <b>Uun</b> 268	432 <b>Uun</b> 268	433 <b>Uun</b> 268	434 <b>U</b>

## Lab 1: TITRATION METHOD

### AIM:

To Titrate Sodium Hydroxide with Hydrochloric acid

### APPARATUS:

Burette (50cm<sup>3</sup>), pipette (10cm<sup>3</sup>), three conical flasks (250cm<sup>3</sup>), two beakers (250cm<sup>3</sup>), funnel, wash bottle, retort stand, boss and clamp, pipette filler.

### MATERIALS:

0.1 mol dm<sup>-3</sup> hydrochloric acid, 0.1 mol dm<sup>-3</sup> sodium hydroxide, phenolphthalein indicator

### PROCEDURE:

- 1) Rinse the burette twice with a little hydrochloric acid and fill up to just above the zero mark.
- 2) Run out acid to remove any air trapped in the jet and refill to above the zero mark.
- 3) Open the tap of the burette and allow acid to run out until the bottom of the meniscus is just on the zero mark when viewed at eye level.
- 4) Remove any hanging drop from the jet by touching against the inside of the beaker containing the acid.
- 5) Draw up sodium hydroxide solution using the pipette to above the graduation mark and allow the solution to run into the sink. This is the rinsing of the pipette.
- 6) Draw up sodium hydroxide solution using the pipette to above the graduation mark and allow the meniscus to fall until the bottom of the meniscus is resting on the graduation mark when viewed at eye level.
- 7) Run the sodium hydroxide solution from the pipette into a conical flask. After all the liquid has run out from the pipette touch the tip of the pipette against the bottom of the flask and withdraw after 15 sec.
- 8) Repeat steps 6 and 7 using the other conical flask.
- 9) Place two drops of phenolphthalein indicator into each flask.
- 10) Perform titrations on the contents of each flask against the acid from the burette.
- 11) Record your results in a table and use this to find the end-point of the titration.

#### Burette Readings

	Rough	Accurate		
		1	2	3
Final burette reading /cm <sup>3</sup>				
Initial burette reading/cm <sup>3</sup>				
Volume used/cm <sup>3</sup>				

### *TREATMENT OF RESULTS:*

- 1) Write a balanced equation for the reaction in this experiment
- 2) Calculate the mean titre volume
- 3) Calculate the number of moles of HCl in the volume used in the titration using your mean titre volume.
- 4) Using your answer to 3 and the mole ratio of the equation, calculate the number of moles of NaOH in the 10.0 cm<sup>3</sup> pipette.
- 5) Calculate the number of moles of NaOH present in 1dm<sup>3</sup> of solution (concentration in moldm<sup>-3</sup>). How close is it to 0.1M?
- 6) Calculate the concentration of NaOH in gdm<sup>-3</sup>

## Lab 2: Determination of Concentration of an Unknown Solution by Titration Method

### AIM:

To determine the concentration of an unknown solution of sodium hydroxide

### APPARATUS:

Burette (50cm<sup>3</sup>), pipette (10cm<sup>3</sup>), three conical flasks (250cm<sup>3</sup>), three beakers (250cm<sup>3</sup>), funnel, wash bottle, retort stand, boss and clamp, pipette filler.

### MATERIALS:

0.1 mol dm<sup>-3</sup> hydrochloric acid, unknown sodium hydroxide solution, phenolphthalein indicator

### PROCEDURE:

- 1) Rinse the burette twice with a little hydrochloric acid and fill up to just above the zero mark.
- 2) Run out acid to remove any air trapped in the jet and refill to above the zero mark.
- 3) Open the tap of the burette and allow acid to run out until the bottom of the meniscus is just on the zero mark when viewed at eye level.
- 4) Remove any hanging drop from the jet by touching against the inside of the beaker containing the acid.
- 5) Draw up sodium hydroxide solution using the pipette to above the graduation mark and allow the solution to run into the waste beaker. This is the rinsing of the pipette.
- 6) Draw up sodium hydroxide solution using the pipette to above the graduation mark and allow the meniscus to fall until the bottom of the meniscus is resting on the graduation mark when viewed at eye level.
- 7) Run the sodium hydroxide solution from the pipette into a conical flask. After all the liquid has run out from the pipette touch the tip of the pipette against the bottom of the flask and withdraw after 15 sec.
- 8) Repeat steps 6 and 7 using the other conical flask.
- 9) Place two drops of phenolphthalein indicator into each flask.
- 10) Perform titrations on each flask against the acid from the burette.
- 11) Record your results in a table and use this to find the end-point of the titration.

#### Burette Readings

	Rough	Accurate		
		1	2	3
Final burette reading /cm <sup>3</sup>				
Initial burette reading/cm <sup>3</sup>				
Volume used/cm <sup>3</sup>				

### *TREATMENT OF RESULTS:*

- 1) Write a balanced equation for the reaction in this experiment
- 2) Calculate the mean titre volume
- 3) Calculate the number of moles of HCl in the volume used in the titration using your mean titre volume.
- 4) Using your answer to 3 and the mole ratio of the equation, calculate the number of moles of NaOH in the 10.0 cm<sup>3</sup> pipette.
- 5) Calculate the number of moles of NaOH present in 1dm<sup>3</sup> of solution (concentration in moldm<sup>-3</sup>)
- 6) Calculate the concentration of NaOH in gdm<sup>-3</sup>

## Lab 3: Redox Titration

### AIM

To Determine the Number of Molecules of Water of Crystallization in ammonium iron (II) sulphate crystals

### APPARATUS


Burette (50cm<sup>3</sup>), pipette (25cm<sup>3</sup>), three conical flasks (250cm<sup>3</sup>), two Beakers (250cm<sup>3</sup>), funnel, wash bottle, retort stand, boss and clamp, pipette filler.

### MATERIAL

0.020 mol dm<sup>-3</sup> potassium manganate(VII) solution, 39.2 g dm<sup>-3</sup> ammonium iron(II) sulphate, 1 mol dm<sup>-3</sup> sulphuric acid.

### THEORY

Ammonium iron (II) sulphate is a double salt with the formula (NH<sub>4</sub>)<sub>2</sub>Fe (SO<sub>4</sub>)<sub>2</sub>.nH<sub>2</sub>O, where n is the number of molecules of water of crystallization. Hence the properties of its solution will be those of its component ions. One of these ions is Fe<sup>2+</sup> and this is the actual reducing agent. The potassium manganate (VII) is the oxidizing agent. The reaction between the MnO<sub>4</sub><sup>-</sup> ion of potassium manganate (VII) and Fe<sup>2+</sup> of the ammonium iron (II) sulphate is given by the redox equation:



From this equation it can be seen that 1 mol of MnO<sub>4</sub><sup>-</sup> reacts with 5 mol of Fe<sup>2+</sup>. Using this equation and the results of the titration, the concentration of the Fe<sup>2+</sup> ions can be determined. From this result, and the information provided, the number of molecules of water of crystallization in the double salt can also be calculated.

### PROCEDURE:

- 1) Rinse the burette twice with a little potassium manganate (VII) solution and fill above the zero mark.
- 2) Run out potassium manganate (VII) to remove any air trapped in the jet and refill to above the zero mark.
- 3) Open the tap of the burette and allow potassium manganate (VII) to run out until the top of the meniscus is just on the zero mark when viewed at eye level. - With manganate (VII) titrations you read the top of the meniscus since the bottom of the meniscus is not clearly visible.
- 4) Remove any hanging drop from the jet by touching against the inside of the beaker containing the potassium manganate (VII).
- 5) Draw up the ammonium iron (II) sulphate solution using the pipette to above the graduation mark and allow the solution to run into the waste beaker. This is the rinsing of the pipette.

- 6) Draw up the ammonium iron (II) sulphate solution using the pipette to above the graduation mark and allow the meniscus to fall until the bottom of the meniscus is resting on the graduation mark when viewed at eye level.
- 7) Run the ammonium iron (II) sulphate solution from the pipette into a conical flask. After all the liquid has run out from the pipette touch the tip of the pipette against the bottom of the flask and withdraw after 15 sec.
- 8) Repeat steps 6 and 7 using the other conical flasks.
- 9) Perform titrations on the contents of each flask against the potassium manganate (VII) from the burette until the first permanent trace of pink is obtained.
- 10) Record your results in a table and repeat the titration until you obtain consistent results.

### RESULTS:

Burette

	Rough	Accurate		
		1	2	3
Final burette reading /cm <sup>3</sup>				
Initial burette reading/cm <sup>3</sup>				
Volume used/cm <sup>3</sup>				

### TREATMENT OF RESULTS:

- 1) Calculate the number of moles of  $\text{MnO}_4^-$  in the volume used in the titration using your mean titre volume.
- 2) Using your answer to 1 and the mole ratio of the equation, calculate the number of moles of  $\text{Fe}^{2+}$  in the 25.0 cm<sup>3</sup> pipette.
- 3) Use your answer to 2 to calculate the concentration of  $\text{Fe}^{2+}$  in mol dm<sup>-3</sup>.
- 4) Divide 39.2 g dm<sup>-3</sup> by your answer to 3 to get the relative molecular mass of the double salt.
- 5)  $284 + 18n = \text{answer from 4}$ . Use this to calculate n.
- 6) Your answer for n to the nearest whole number is the number of molecules of water of crystallization in ammonium iron (II) sulphate.

# **APPENDIX A**

## **LAB REPORT SHEET FORMAT**

**NAME:**

**PARTNER(S) NAME:**

**DATE:**

**TITLE:**

**AIM:** brief description of the purpose of the lab.

**APPARATUS:** list of all the equipment and chemicals that were used in the exercise.

**METHOD:** written in past perfect tense, this is a numbered sequence of the steps that were followed to carry out the exercise.

For example: The addition of iodine to a test tube would be described in the following way:

“The iodine **was added** to the test tube...”

NOT

“They added iodine to the test tube...”

### **PRECAUTIONS FOR ACCURATE RESULTS AND SAFETY:**

- List of two important precautions which must be followed when using the equipment for the lab exercise.
- For the procedure followed in the virtual lab, suggest two safety measures which must have been observed.

**RESULTS:** present the results of the experiment (the things you saw or heard). Draw a table wherever possible (title should be written in caps and underlined **above** the table).

**CONCLUSION:** in relation to the aim of the experiment, the deductions made based on the results obtained.

**ANSWERS TO QUESTIONS:** answer any post lab questions.

## **APPENDIX B**

### **EXAMPLE OF A COMPLETED LAB REPORT SHEET**

**NAME:** Jane Joelly

**PARTNER'S NAME:** Richard Rally

**DATE:** 20.05.2016

**TITLE:** Osmosis

**AIM:** To observe the process of osmosis.

**APPARATUS:** 3 cups, vinegar, corn syrup, water, egg, stop-clock.

#### **METHOD:**

1. The egg was soaked in vinegar for two days.
2. The egg shell was then removed leaving the soft membrane still holding the yolk and albumen.
3. The egg was then placed in a cup containing enough corn syrup to cover most of it and left for an hour.
4. Any difference in size and appearance was recorded.
5. The egg was then placed in another cup containing just enough fresh water to cover most of it for another two hours.
6. Any difference in size and appearance was recorded.

#### **PRECAUTIONS FOR ACCURATE RESULTS AND SAFETY:**

Accuracy:

- The stop clock was tested before use to ensure that it was working properly.
- The batteries in the stopclock were changed to ensure that they did not run out during the experiment.

Safety:

- Gloves should have been worn to protect the hands from vinegar.
- Protective glasses were worn to protect the eyes from the vinegar.

#### **RESULTS:**

After 1 hour in the corn syrup, the egg was smaller than it was at first. There were also wrinkles on the surface.

After 10 minutes in the fresh water, the wrinkles disappeared. The longer the egg stayed in the fresh water, the larger it swelled.

## **CONCLUSION:**

In this experiment, osmosis was observed. Osmosis is the movement of water particles from a dilute solution to a more concentrated solution across a semi-permeable membrane. When the egg shell was removed, the membrane left behind was semi-permeable. The albumen was more dilute than the corn syrup, so water moved from the albumen across the semi-permeable membrane and into the syrup. This caused the egg to reduce in size and wrinkle. The albumen was more concentrated than the fresh water however, so water particles moved across the semi-permeable membrane into the egg, making it swell.

## **ANSWERS TO QUESTIONS:**

1. Would the same observations be made with an egg with the shell intact?

No, because the egg shell is not semi-permeable.

2. If the experiment is carried out using the same method but replacing the corn syrup with salt water what observations would be made?

The same observations would be made once the concentration of the salt solution was higher than that of the albumen.

## APPENDIX C

### LAB RUBRIC

This is the Mark Scheme which will be used to score your Lab Reports

SECTION	COMPETENCY	PERFORMANCE RATING			
		Excellent 4	Average 3	Fair 2	Poor 1
<b>Aim</b>	Effective Communication	The aim is clearly stated.	The aim is slightly unclear.	The aim is quite vague.	The aim is not stated.
<b>Apparatus</b>	Knowledge	Apparatus list is complete. Student gives the correct name of all equipment used.	Apparatus list is complete. Student gives the correct name of most equipment used.	Apparatus list is incomplete. Student gives the correct name of most the equipment they listed.	Apparatus list is incomplete. Student gives the correct name of only a few pieces of the equipment they listed.
<b>Method</b>	Effective Communication	Method is correct, logical, complete and written in appropriate tense.	Method is incorrect, or illogical, or incomplete or written in inappropriate tense.	Method has two - three deficiencies.	Method is incorrect and illogical and incomplete and written in inappropriate tense.
	Technical expertise	Students suggest at least two precautions and two safety measures for each lab which are relevant and accurate.	Students suggest at least one precaution and one safety measure for each lab is relevant and accurate.	Students suggest at least one precaution or one safety measure for each lab which is relevant and accurate.	Students do not suggest any precaution or safety measure which is relevant and accurate.

<b>Results</b>	Effective Communication	Each observation is accurately and completely included.	Most observations are accurately and completely included.	Few observations are accurately and completely included.	No observations are accurately and completely included.
	Technological Expertise	All results are recorded in appropriate format. All calculations, tables and graphs are done correctly.	Most results are recorded in appropriate format or most calculations, tables and graphs are done correctly.	Most results are not recorded in appropriate format or most calculations, tables and graphs are not done correctly.	Most results are not recorded in appropriate format and most calculations, tables and graphs are not done correctly.
<b>Conclusion</b>	Analytical Ability	The chemical principles underlying each observation are fully explained and relevant conclusions are made.	The chemical principles underlying most observations are fully explained and relevant conclusions are made.	The chemical principles underlying few observations are fully explained and relevant conclusions are made.	No observations are fully explained and irrelevant conclusions are made.
<b>Post lab questions/ Post Lab Worksheet</b>	Analytical Ability	All questions answered correctly.	Most questions answered correctly.	Few questions answered correctly.	No questions answered correctly.

## APPENDIX D

### LAB ASSESSMENT SHEET

The marks for each lab you submit will be recorded on one of these sheets by the Lab Lecturer.

<b>i. Student Name:</b>						
<b>ii. Course:</b>						
<b>iii. CRN:</b>						
<b>iv. Lab #:</b>						
<b>v. Date:</b>						
Section	Competency	SCORE				COMMENTS
<b>Aim</b>	Eff. Comm. 1	4	3	2	1	
	<b>TOTAL</b>					
<b>Apparatus</b>	Knowledge	4	3	2	1	
	<b>TOTAL</b>					
<b>Method</b>	Eff. Comm.	4	3	2	1	
	Tech. Exp.	4	3	2	1	
	<b>TOTAL</b>					
<b>Results</b>	Tech. Exp.	4	3	2	1	
	Eff. Comm.	4	3	2	1	
	<b>TOTAL</b>					
<b>Conclusion</b>	Analytical Ability	4	3	2	1	
	<b>TOTAL</b>					
<b>Post Lab questions</b>	Analytical Ability	4	3	2	1	
	<b>TOTAL</b>					
<b>TOTAL LAB MARK /32:</b>						
<b>TOTAL LAB MARK (%):</b>						
<b>LAB LECTURER:</b>						