



**COLLEGE OF SCIENCE, TECHNOLOGY & APPLIED ARTS**

**OF TRINIDAD & TOBAGO**

**School of Nursing Health and Environmental Sciences  
Department of Natural & Life Sciences**


---

## **CHEM 090 – Introduction to Concepts in Chemistry**

---

**Student Manual**

Edited: 29/02/2016



## A Word of Welcome

The School of Nursing Health and Environmental Sciences and the Department of Natural & Life Sciences welcome you to the practical component of the course. The experiments that form part of this course have been selected to:

- Reinforce your knowledge and understanding of the theoretical content of your first year Chemistry course.
- Develop skills such as observing, recording, measuring and reporting of results.
- Teach you special techniques, such as weighing, pipetting, titrating, filtration etc, so that you will become familiar with the basics of laboratory work.



## TABLE OF CONTENTS

A WORD OF WELCOME .....	2
TABLE OF CONTENTS .....	3
INTRODUCTION .....	4
SAFETY IN THE CHEMISTRY LABORATORY.....	7
WASTE DISPOSAL.....	13
MANAGING YOUR TIME DURING THE PRACTICAL.....	14
GETTING STARTED.....	15
LAB 1A: INTRODUCTION TO CHEMISTRY PRACTICALS .....	18
LAB 1B: DETERMINATION OF MELTING POINT .....	20
LAB 2: SEPARATION OF MIXTURES.....	22
LAB 3: EFFECT OF TEMPERATURE ON THE SOLUBILITY OF A SALT .....	26

## Introduction

### General

Always read and study the relevant experiment before coming to the laboratory. Check out appropriate references to help you understand the procedures you will be carrying out.

### **Attendance at all laboratory sessions is compulsory.**

Learning lab techniques must come from *experiencing* lab techniques. Experience is gained from carrying out experiments. You will need to plan ahead for your lab session, and think about what you are doing (and why) while you are in lab.

### **Use past experiences to develop new skills.**

You will gain new skills with new equipment and learn more about the chemistry theories and concepts you hear about in lectures, but the skills that we expect you to develop the most during these experiments are general scientific skills that should be transferable to any other scientific discipline and to your daily life. These skills developed in the laboratory are fundamental to the world of work because industries depend on the laboratory to maintain a level of quality, which enhances the company's reputation, as well as for product development.

The lab is a part of the learning process.

### **LABORATORY SCHEDULE:**

The Lab schedules are posted at the start of the semester on notice boards, lab tech stations and on the lab doors. Please refer to them to find out the dates and times of your lab sessions.

Once the semester begins, you will have to balance your laboratory study time with all of the other demands that you have. Most of the experiments have been scheduled so that the average student can finish everything (including the report) required in the laboratory period.

**You are expected to attend all scheduled laboratory periods and make-up labs will not be scheduled.**

# SONHES Multi-Science Laboratory Rules

---

## Before Entering the Lab

**Proper attire is required for entry into the lab\***

- Lab Coat
  - Long sleeve (wrist length)
  - Knee Length
  - Drill cotton
- Long pants or long skirt (ankle length)
- Shoes that fully cover your feet
- Long hair tied back
- No dangling jewelry or accessories

## Punctuality

- Students are not allowed entry into lab 15 minutes **after** the start of the lab session.


## During the Lab Session

- No eating, No drinking, No chewing of gum, No smoking
- Mobile phones and electronic devices (laptops, tablets etc.) are not to be used in the lab.
- Students are expected to follow the lab rules and to follow the guidance of the lab technicians while in the lab.
- Students are expected to have their own lab manuals at each lab session. No new manuals will be issued to students who forgot to bring theirs.

## Leaving the Lab

- Students are required to wash all glassware, return items to student cupboards as well as dispose of gloves, tissue, scrap paper etc and wipe countertops before leaving the lab.
- Lab coats are to be taken off when leaving the lab (including trips to the bathroom).


**\*LAB LECTURERS OR LAB TECHNICIANS MAY PREVENT ENTRY OR ASK STUDENTS TO LEAVE THE LAB IN THE EVENT OF STUDENT NON-COMPLIANCE TO THE LABORATORY RULES.**



- ✓ Wrist length sleeves for coverage
- ✓ Protects clothes
- ✓ Able to button at front
- ✓ Thick material for protection from spills
- ✓ Durable, long lasting protection for future labs
- ✓ Feel like a scientist!!



- ✗ Minimal barrier against chemicals and microbes
- ✗ Skin exposed
- ✗ Thin material allows chemicals to permeate and reach skin
- ✗ Hoodies can transmit contaminants to friends & family




- ✓ Feet fully covered!
- ✓ Comfortable
- ✓ Ankle length Pants
- ✓ Protection from spills and broken glassware



- ✗ Feet exposed
- ✗ Limited protection from spills
- ✗ Possibility of cuts from broken glassware
- ✗ Possible exposure to pathogens

## SAFETY IN THE CHEMISTRY LABORATORY

The laboratory is a hazardous environment in which to work. The hazards are often unavoidable since scientists regularly have to use hazardous materials and equipment. However, with sensible precautions the laboratory is probably no more dangerous than your home. You **MUST AT ALL TIMES OBSERVE ALL RULES** stated in this manual.



## Proper Lab wear

Gloves, lab coat, proper lab shoes

Barrier to prevent contact of contaminants with skin

## Be Punctual

Check the schedule regularly, arrive to lab on time

Be present for all safety information regarding the lab

## Good Hygiene

Wash hands when leaving the lab, swab counter tops after lab

Prevents cross contamination to surfaces outside of lab

## Proper use of Equipment

Check electric cords before use

Always use equipment for its intended purpose



## Essential Rules for Laboratory Safety

The essential rules for laboratory safety can be expressed under two simple headings:

### **ALWAYS AND NEVER**

#### **ALWAYS**

- ♦ Familiarize yourself with the laboratory safety procedures.
- ♦ Dress appropriately.
- ♦ Wash your hands before leaving the laboratory.
- ♦ Read the instructions for the laboratory protocol **carefully** before starting any experiment.
- ♦ Handle all chemicals with great care.
- ♦ Keep your working area tidy.
- ♦ Immediately inform your instructor about spills.

#### **NEVER**

- ♦ Eat or drink in the laboratory. CHEWING OF GUM IS STRICTLY PROHIBITED!
- ♦ Smoke in the laboratory
- ♦ Inhale, taste or sniff chemicals.
- ♦ Fool around or distract neighbours.
- ♦ Carry out unauthorized experiments.

### **Laboratory Safety Procedures**

Your laboratory will have certain procedures, which you must be familiar with. Make sure you know where all the exits from the laboratory are, in the event of an evacuation because of fire or other incident. Locate the SHOWERS and eyewash stations and know the type the fire extinguishers that are in the laboratory and how to operate them.

### **Dress Code**

- Wear a Lab coat.
- Wear shoes that cover your entire feet (sneakers) Avoid open footwear & high heels. So sandals, slipper and ballet flats are unacceptable.

- Wear approved safety goggles or safety glasses when in the lab.
- Avoid floppy garments; avoid things that dangle. These get tangled up in equipment or glassware and cause accidents.
- Avoid long, loose hair styles for the same reason. When Bunsen burners are in use, long hair sometimes catches on fire.

## Laboratory Techniques

### Handling of Acids or Bases

- In addition to the dress code above, personal protective equipment, such as Vinyl" gloves or "nitrile" GLOVES SHOULD be worn.
- For even more dangerous liquids, thick "nitrile" gloves should be used. These are very safe.
- **DO NOT add WATER to CONCENTRATED ACID.** The heat generated will cause splattering. If necessary to prepare certain solutions, **DO** add ACID to WATER (instead of the reverse order of addition). The heat generated will be less, but splattering still may occur. A good practice in all lab operations is to keep things at arm's length.

### Become acquainted with handling laboratory safety equipment e.g.

- Fire Extinguishers
- Fire Blanket
- Eye-wash Fountain
- First-Aid Kit

### Become acquainted with using of the FUME HOOD

- Use the HOOD for reactions that give off vapours, especially smelly vapours.
- The draft of the HOOD will sweep away vapours so that the lab itself maintains reasonable air quality.

### Locate the Safety Shower and understand how it works

- Shower should be used for **dire EMERGENCY** only!
- If you are **ON FIRE**, or suffer a **massive spill** of a **dangerous chemical**, and need to get it off rapidly.

**Position yourself under the safety shower and pull the handle—a deluge of water will result.**

**Immediately clean up Broken Glass**

- Sweep it up right away
- Place the broken glass in a "SHARP'S CONTAINER."
- This is a thick walled carton that will be sealed and discarded as such.


**Student Use of the Laboratory**

- No UNAUTHORIZED experiments! These are terrible dangers in unskilled hands
- No EATING or DRINKING in the LAB. A good practice is to assume everything in the lab is toxic.
- DO NOT drink from lab equipment. NOT EVEN distilled water.
- When using pipettes, **DO NOT** suck the liquid into the pipette directly by mouth.
- Frequent mistakes, even by skilled lab workers, lead to the liquid reaching the mouth.
- Do NOT smell the flask directly.

**Housekeeping:**


- **STUDENTS MUST CLEAN UP AFTER THEMSELVES.**
- Paper and trash must not be left in the room, thrown into drawers, sinks, on the floor, etc.
- When the equipment boxes are on carts, students who take a box from the cart to do an experiment must put everything back into the box neatly and put the box back on the cart.

If you absolutely must test the odour, carefully waft the vapour from the flask toward your nose with your hand, keeping the flask quite distant from your face.



***Correct way of wafting odors.***

Make sure that reaction tubes, e.g. test tubes, are not directed toward yourself or other persons. The chemicals may splatter out the tube.




## WASTE DISPOSAL

We all have a responsibility towards a clean and safe environment. The Chemistry Department's policy dealing with the disposal of laboratory waste is in accordance with the Occupational Health and Safety Act, which is a law of our country. We do not expect you to understand every aspect of this law, but as responsible citizen and laboratory user you must be aware of the correct procedures for disposal of the different classes of laboratory waste. You will be given specific instructions on how to dispose of your chemical wastes during each pre-lab briefing. Here are also some general rules that you must take note of:

1. Place broken glass into the specially marked bins.
2. Drain all harmless chemicals in solution into the sink followed by plenty of running water.
3. Pour toxic inorganic waste solutions into the specially labelled container in the fume cupboard.
4. Discard paper and any other solid waste into the bin.
5. Ensure that matches are extinguished before disposing of them in this way.
6. Shut off all gas and water lines when not in use.

## MANAGING YOUR TIME DURING THE PRACTICAL

1. It is recommended that you arrange your workbench according to the diagram below. Keeping your workspace organized and free of clutter will save you time and frustration, and contributes greatly to laboratory safety.



2. When collecting chemicals:
  - Choose the size of your container according to the volume of chemical that you will be collecting, for instance: for 100 cm<sup>3</sup> of solution, use a 250 cm<sup>3</sup> beaker and for 10 cm<sup>3</sup> use a test tube.
  - Before collecting chemicals and/or solutions from the dispensary, mark each container with the name or formula of the chemical to be collected. This will prevent any mix-ups later on.
3. Use a small notebook for jotting down masses, measurements and observations. Writing these bits of information on slips of scrap paper that can easily be lost is not only unprofessional, but also risky because it means that your whole afternoon's practical work can be wasted because of lost data.

## GETTING STARTED

For each laboratory period, including the first, you will be required to:

- ◆ Study the entire experimental outline in the **Laboratory Manual**. It includes specific directions concerning laboratory philosophy and protocols.
- ◆ Prepare your **Laboratory Notebook**. Include a summary and a procedural outline for your experiment (see "Rules for Keeping Your Laboratory Notebook" and Instruction in the experiment outline).
- ◆ Bring your safety goggles or use the safety goggles provided. These offer adequate protection against accidental splashing of corrosive chemicals.
- ◆ Come to lab dressed appropriately

### A GOOD PRACTICE:


Read the experimental procedure ahead of lab.

Avoid horseplay. In a laboratory setting, horseplay, even if good-natured, is absolutely unacceptable.

- **No pushing!**
- **No shoving!**

At the end of the lab period: Exit the lab in an orderly manner. Again: no running, no pushing, AND NO shoving.

Figure 1: Apparatus that you will be using in the laboratory





Atomic number																			
Symbol																			
Atomic weight																			
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18		
1 <b>H</b> 1.008	4 <b>Be</b> 9.012	9 <b>Li</b> 6.941	12 <b>Mg</b> 24.31	19 <b>K</b> 39.10	20 <b>Ca</b> 40.08	23 <b>V</b> 50.94	24 <b>Cr</b> 52.00	25 <b>Mn</b> 54.94	26 <b>Fe</b> 55.85	27 <b>Co</b> 58.93	28 <b>Ni</b> 58.69	29 <b>Cu</b> 63.55	30 <b>Zn</b> 65.39	31 <b>Ga</b> 69.72	32 <b>Ge</b> 72.61	33 <b>As</b> 74.92	34 <b>Se</b> 78.96	35 <b>Br</b> 79.90	36 <b>Kr</b> 83.80
37 <b>Rb</b> 85.47	38 <b>Sr</b> 87.62	39 <b>Y</b> 88.91	40 <b>Zr</b> 91.22	41 <b>Nb</b> 92.91	42 <b>Mo</b> 95.94	43 <b>Tc</b> 98.91	44 <b>Ru</b> 101.1	45 <b>Rh</b> 102.9	46 <b>Pd</b> 106.4	47 <b>Ag</b> 107.9	48 <b>Cd</b> 112.4	49 <b>In</b> 114.8	50 <b>Sn</b> 118.7	51 <b>Sb</b> 121.8	52 <b>Te</b> 127.6	53 <b>I</b> 126.9	54 <b>Xe</b> 131.3	55 <b>Cs</b> 132.9	56 <b>Ba</b> 137.3
87 <b>Fr</b> 223.0	88 <b>Ra</b> 226.0	89 <b>Ac</b> 227.0	90 <b>Th</b> 232.0	91 <b>Pa</b> 231.0	92 <b>U</b> 238.0	93 <b>Np</b> 237.0	94 <b>Pu</b> 244.1	95 <b>Am</b> 243.1	96 <b>Cm</b> 247.1	97 <b>Bk</b> 247.1	98 <b>Cf</b> 251.1	99 <b>Es</b> 252.0	100 <b>Fm</b> 257.1	101 <b>Md</b> 258.1	102 <b>No</b> 259.1	103 <b>Lr</b> 262.1	104 <b>Hf</b> 178.5	105 <b>Ta</b> 180.9	106 <b>W</b> 183.8
103 <b>Lr</b> 262.1	104 <b>Hf</b> 178.5	105 <b>Ta</b> 180.9	106 <b>W</b> 183.8	107 <b>Re</b> 186.2	108 <b>Os</b> 190.2	109 <b>Ir</b> 192.2	110 <b>Pt</b> 195.1	111 <b>Au</b> 197.0	112 <b>Hg</b> 200.6	113 <b>Tl</b> 204.4	114 <b>Pb</b> 207.2	115 <b>Bi</b> 209.0	116 <b>Po</b> 209.0	117 <b>At</b> 210.0	118 <b>Rn</b> 222.0	119 <b>Uuo</b> 293	120 <b>Uuh</b> 289	121 <b>Uus</b> 293	122 <b>Uut</b> 289
123 <b>Uut</b> 289	124 <b>Uuq</b> 289	125 <b>Uup</b> 289	126 <b>Uuh</b> 289	127 <b>Uus</b> 293	128 <b>Uuo</b> 293	129 <b>Uuh</b> 289	130 <b>Uus</b> 293	131 <b>Uut</b> 289	132 <b>Uuq</b> 289	133 <b>Uup</b> 289	134 <b>Uuh</b> 289	135 <b>Uus</b> 293	136 <b>Uuo</b> 293	137 <b>Uuh</b> 289	138 <b>Uus</b> 293	139 <b>Uuo</b> 293	140 <b>Uuh</b> 289	141 <b>Uus</b> 293	142 <b>Uuo</b> 293
143 <b>Uuh</b> 289	144 <b>Uus</b> 293	145 <b>Uuo</b> 293	146 <b>Uuh</b> 289	147 <b>Uus</b> 293	148 <b>Uuo</b> 293	149 <b>Uuh</b> 289	150 <b>Uus</b> 293	151 <b>Uuo</b> 293	152 <b>Uuh</b> 289	153 <b>Uus</b> 293	154 <b>Uuo</b> 293	155 <b>Uuh</b> 289	156 <b>Uus</b> 293	157 <b>Uuo</b> 293	158 <b>Uuh</b> 289	159 <b>Uus</b> 293	160 <b>Uuo</b> 293	161 <b>Uuh</b> 289	162 <b>Uus</b> 293
163 <b>Uuh</b> 289	164 <b>Uus</b> 293	165 <b>Uuo</b> 293	166 <b>Uuh</b> 289	167 <b>Uus</b> 293	168 <b>Uuo</b> 293	169 <b>Uuh</b> 289	170 <b>Uus</b> 293	171 <b>Uuo</b> 293	172 <b>Uuh</b> 289	173 <b>Uus</b> 293	174 <b>Uuo</b> 293	175 <b>Uuh</b> 289	176 <b>Uus</b> 293	177 <b>Uuo</b> 293	178 <b>Uuh</b> 289	179 <b>Uus</b> 293	180 <b>Uuo</b> 293	181 <b>Uuh</b> 289	182 <b>Uus</b> 293
183 <b>Uuh</b> 289	184 <b>Uus</b> 293	185 <b>Uuo</b> 293	186 <b>Uuh</b> 289	187 <b>Uus</b> 293	188 <b>Uuo</b> 293	189 <b>Uuh</b> 289	190 <b>Uus</b> 293	191 <b>Uuo</b> 293	192 <b>Uuh</b> 289	193 <b>Uus</b> 293	194 <b>Uuo</b> 293	195 <b>Uuh</b> 289	196 <b>Uus</b> 293	197 <b>Uuo</b> 293	198 <b>Uuh</b> 289	199 <b>Uus</b> 293	200 <b>Uuo</b> 293	201 <b>Uuh</b> 289	202 <b>Uus</b> 293
203 <b>Uuh</b> 289	204 <b>Uus</b> 293	205 <b>Uuo</b> 293	206 <b>Uuh</b> 289	207 <b>Uus</b> 293	208 <b>Uuo</b> 293	209 <b>Uuh</b> 289	210 <b>Uus</b> 293	211 <b>Uuo</b> 293	212 <b>Uuh</b> 289	213 <b>Uus</b> 293	214 <b>Uuo</b> 293	215 <b>Uuh</b> 289	216 <b>Uus</b> 293	217 <b>Uuo</b> 293	218 <b>Uuh</b> 289	219 <b>Uus</b> 293	220 <b>Uuo</b> 293	221 <b>Uuh</b> 289	222 <b>Uus</b> 293
223 <b>Uuh</b> 289	224 <b>Uus</b> 293	225 <b>Uuo</b> 293	226 <b>Uuh</b> 289	227 <b>Uus</b> 293	228 <b>Uuo</b> 293	229 <b>Uuh</b> 289	230 <b>Uus</b> 293	231 <b>Uuo</b> 293	232 <b>Uuh</b> 289	233 <b>Uus</b> 293	234 <b>Uuo</b> 293	235 <b>Uuh</b> 289	236 <b>Uus</b> 293	237 <b>Uuo</b> 293	238 <b>Uuh</b> 289	239 <b>Uus</b> 293	240 <b>Uuo</b> 293	241 <b>Uuh</b> 289	242 <b>Uus</b> 293
243 <b>Uuh</b> 289	244 <b>Uus</b> 293	245 <b>Uuo</b> 293	246 <b>Uuh</b> 289	247 <b>Uus</b> 293	248 <b>Uuo</b> 293	249 <b>Uuh</b> 289	250 <b>Uus</b> 293	251 <b>Uuo</b> 293	252 <b>Uuh</b> 289	253 <b>Uus</b> 293	254 <b>Uuo</b> 293	255 <b>Uuh</b> 289	256 <b>Uus</b> 293	257 <b>Uuo</b> 293	258 <b>Uuh</b> 289	259 <b>Uus</b> 293	260 <b>Uuo</b> 293	261 <b>Uuh</b> 289	262 <b>Uus</b> 293
263 <b>Uuh</b> 289	264 <b>Uus</b> 293	265 <b>Uuo</b> 293	266 <b>Uuh</b> 289	267 <b>Uus</b> 293	268 <b>Uuo</b> 293	269 <b>Uuh</b> 289	270 <b>Uus</b> 293	271 <b>Uuo</b> 293	272 <b>Uuh</b> 289	273 <b>Uus</b> 293	274 <b>Uuo</b> 293	275 <b>Uuh</b> 289	276 <b>Uus</b> 293	277 <b>Uuo</b> 293	278 <b>Uuh</b> 289	279 <b>Uus</b> 293	280 <b>Uuo</b> 293	281 <b>Uuh</b> 289	282 <b>Uus</b> 293
283 <b>Uuh</b> 289	284 <b>Uus</b> 293	285 <b>Uuo</b> 293	286 <b>Uuh</b> 289	287 <b>Uus</b> 293	288 <b>Uuo</b> 293	289 <b>Uuh</b> 289	290 <b>Uus</b> 293	291 <b>Uuo</b> 293	292 <b>Uuh</b> 289	293 <b>Uus</b> 293	294 <b>Uuo</b> 293	295 <b>Uuh</b> 289	296 <b>Uus</b> 293	297 <b>Uuo</b> 293	298 <b>Uuh</b> 289	299 <b>Uus</b> 293	300 <b>Uuo</b> 293	301 <b>Uuh</b> 289	302 <b>Uus</b> 293
303 <b>Uuh</b> 289	304 <b>Uus</b> 293	305 <b>Uuo</b> 293	306 <b>Uuh</b> 289	307 <b>Uus</b> 293	308 <b>Uuo</b> 293	309 <b>Uuh</b> 289	310 <b>Uus</b> 293	311 <b>Uuo</b> 293	312 <b>Uuh</b> 289	313 <b>Uus</b> 293	314 <b>Uuo</b> 293	315 <b>Uuh</b> 289	316 <b>Uus</b> 293	317 <b>Uuo</b> 293	318 <b>Uuh</b> 289	319 <b>Uus</b> 293	320 <b>Uuo</b> 293	321 <b>Uuh</b> 289	322 <b>Uus</b> 293
323 <b>Uuh</b> 289	324 <b>Uus</b> 293	325 <b>Uuo</b> 293	326 <b>Uuh</b> 289	327 <b>Uus</b> 293	328 <b>Uuo</b> 293	329 <b>Uuh</b> 289	330 <b>Uus</b> 293	331 <b>Uuo</b> 293	332 <b>Uuh</b> 289	333 <b>Uus</b> 293	334 <b>Uuo</b> 293	335 <b>Uuh</b> 289	336 <b>Uus</b> 293	337 <b>Uuo</b> 293	338 <b>Uuh</b> 289	339 <b>Uus</b> 293	340 <b>Uuo</b> 293	341 <b>Uuh</b> 289	342 <b>Uus</b> 293
343 <b>Uuh</b> 289	344 <b>Uus</b> 293	345 <b>Uuo</b> 293	346 <b>Uuh</b> 289	347 <b>Uus</b> 293	348 <b>Uuo</b> 293	349 <b>Uuh</b> 289	350 <b>Uus</b> 293	351 <b>Uuo</b> 293	352 <b>Uuh</b> 289	353 <b>Uus</b> 293	354 <b>Uuo</b> 293	355 <b>Uuh</b> 289	356 <b>Uus</b> 293	357 <b>Uuo</b> 293	358 <b>Uuh</b> 289	359 <b>Uus</b> 293	360 <b>Uuo</b> 293	361 <b>Uuh</b> 289	362 <b>Uus</b> 293
363 <b>Uuh</b> 289	364 <b>Uus</b> 293	365 <b>Uuo</b> 293	366 <b>Uuh</b> 289	367 <b>Uus</b> 293	368 <b>Uuo</b> 293	369 <b>Uuh</b> 289	370 <b>Uus</b> 293	371 <b>Uuo</b> 293	372 <b>Uuh</b> 289	373 <b>Uus</b> 293	374 <b>Uuo</b> 293	375 <b>Uuh</b> 289	376 <b>Uus</b> 293	377 <b>Uuo</b> 293	378 <b>Uuh</b> 289	379 <b>Uus</b> 293	380 <b>Uuo</b> 293	381 <b>Uuh</b> 289	382 <b>Uus</b> 293
383 <b>Uuh</b> 289	384 <b>Uus</b> 293	385 <b>Uuo</b> 293	386 <b>Uuh</b> 289	387 <b>Uus</b> 293	388 <b>Uuo</b> 293	389 <b>Uuh</b> 289	390 <b>Uus</b> 293	391 <b>Uuo</b> 293	392 <b>Uuh</b> 289	393 <b>Uus</b> 293	394 <b>Uuo</b> 293	395 <b>Uuh</b> 289	396 <b>Uus</b> 293	397 <b>Uuo</b> 293	398 <b>Uuh</b> 289	399 <b>Uus</b> 293	400 <b>Uuo</b> 293	401 <b>Uuh</b> 289	402 <b>Uus</b> 293
403 <b>Uuh</b> 289	404 <b>Uus</b> 293	405 <b>Uuo</b> 293	406 <b>Uuh</b> 289	407 <b>Uus</b> 293	408 <b>Uuo</b> 293	409 <b>Uuh</b> 289	410 <b>Uus</b> 293	411 <b>Uuo</b> 293	412 <b>Uuh</b> 289	413 <b>Uus</b> 293	414 <b>Uuo</b> 293	415 <b>Uuh</b> 289	416 <b>Uus</b> 293	417 <b>Uuo</b> 293	418 <b>Uuh</b> 289	419 <b>Uus</b> 293	420 <b>Uuo</b> 293	421 <b>Uuh</b> 289	422 <b>Uus</b> 293
423 <b>Uuh</b> 289	424 <b>Uus</b> 293	425 <b>Uuo</b> 293	426 <b>Uuh</b> 289	427 <b>Uus</b> 293	428 <b>Uuo</b> 293	429 <b>Uuh</b> 289	430 <b>Uus</b> 293	431 <b>Uuo</b> 293	432 <b>Uuh</b> 289	433 <b>Uus</b> 293	434 <b>Uuo</b> 293	435 <b>Uuh</b> 289	436 <b>Uus</b> 293	437 <b>Uuo</b> 293	438 <b>Uuh</b> 289	439 <b>Uus</b> 293	440 <b>Uuo</b> 293	441 <b>Uuh</b> 289	442 <b>Uus</b> 293
443 <b>Uuh</b> 289	444 <b>Uus</b> 293	445 <b>Uuo</b> 293	446 <b>Uuh</b> 289	447 <b>Uus</b> 293	448 <b>Uuo</b> 293	449 <b>Uuh</b> 289	450 <b>Uus</b> 293	451 <b>Uuo</b> 293	452 <b>Uuh</b> 289	453 <b>Uus</b> 293	454 <b>Uuo</b> 293	455 <b>Uuh</b> 289	456 <b>Uus</b> 293	457 <b>Uuo</b> 293	458 <b>Uuh</b> 289	459 <b>Uus</b> 293	460 <b>Uuo</b> 293	461 <b>Uuh</b> 289	462 <b>Uus</b> 293
463 <b>Uuh</b> 289	464 <b>Uus</b> 293	465 <b>Uuo</b> 293	466 <b>Uuh</b> 289	467 <b>Uus</b> 293	468 <b>Uuo</b> 293	469 <b>Uuh</b> 289	470 <b>Uus</b> 293	471 <b>Uuo</b> 293	472 <b>Uuh</b> 289	473 <b>Uus</b> 293	474 <b>Uuo</b> 293	475 <b>Uuh</b> 289	476 <b>Uus</b> 293	477 <b>Uuo</b> 293	478 <b>Uuh</b> 289	479 <b>Uus</b> 293	480 <b>Uuo</b> 293	481 <b>Uuh</b> 289	482 <b>Uus</b> 293
483 <b>Uuh</b> 289	484 <b>Uus</b> 293	485 <b>Uuo</b> 293	486 <b>Uuh</b> 289	487 <b>Uus</b> 293	488 <b>Uuo</b> 293	489 <b>Uuh</b> 289	490 <b>Uus</b> 293	491 <b>Uuo</b> 293	492 <b>Uuh</b> 289	493 <b>Uus</b> 293	494 <b>Uuo</b> 293	495 <b>Uuh</b> 289	496 <b>Uus</b> 293	497 <b>Uuo</b> 293	498 <b>Uuh</b> 289	499 <b>Uus</b> 293	500 <b>Uuo</b> 293	501 <b>Uuh</b> 289	502 <b>Uus</b> 293
503 <b>Uuh</b> 289	504 <b>Uus</b> 293	505 <b>Uuo</b> 293	506 <b>Uuh</b> 289	507 <b>Uus</b> 293	508 <b>Uuo</b> 293	509 <b>Uuh</b> 289	510 <b>Uus</b> 293	511 <b>Uuo</b> 293	512 <b>Uuh</b> 289	513 <b>Uus</b> 293	514 <b>Uuo</b> 293	515 <b>Uuh</b> 289	516 <b>Uus</b> 293	517 <b>Uuo</b> 293	518 <b>Uuh</b> 289	519 <b>Uus</b> 293	520 <b>Uuo</b> 293	521 <b>Uuh</b> 289	522 <b>Uus</b> 293
523 <b>Uuh</b> 289	524 <b>Uus</b> 293	525 <b>Uuo</b> 293	526 <b>Uuh</b> 289	527 <b>Uus</b> 293	528 <b>Uuo</b> 293	529 <b>Uuh</b> 289	530 <b>Uus</b> 293	531 <b>Uuo</b> 293	532 <b>Uuh</b> 289	533 <b>Uus</b> 293	534 <b>Uuo</b> 293	535 <b>Uuh</b> 289	536 <b>Uus</b> 293	537 <b>Uuo</b> 293	538 <b>Uuh</b> 289	539 <b>Uus</b> 293	540 <b>Uuo</b> 293	541 <b>Uuh</b> 289	542 <b>Uus</b> 293
543 <b>Uuh</b> 289	544 <b>Uus</b> 293	545 <b>Uuo</b> 293	546 <b>Uuh</b> 289	547 <b>Uus</b> 293	548 <b>Uuo</b> 293	549 <b>Uuh</b> 289	550 <b>Uus</b> 293	551 <b>Uuo</b> 293	552 <b>Uuh</b> 289	553 <b>Uus</b> 293	554 <b>Uuo</b> 293	555 <b>Uuh</b> 289	556 <b>Uus</b> 293	557 <b>Uuo</b> 293	558 <b>Uuh</b> 289	559 <b>Uus</b> 293	560 <b>Uuo</b> 293	561 <b>Uuh</b> 289	562 <b>Uus</b> 293
563 <b>Uuh</b> 289	564 <b>Uus</b> 293	565 <b>Uuo</b> 293	566 <b>Uuh</b> 289	567 <b>Uus</b> 293	568 <b>Uuo</b> 293	569 <b>Uuh</b> 289	570 <b>Uus</b> 293	571 <b>Uuo</b> 293	572 <b>Uuh</b> 289	573 <b>Uus</b> 293	574 <b>Uuo</b> 293	575 <b>Uuh</b> 289	576 <b>Uus</b> 293	577 <b>Uuo</b> 293	578 <b>Uuh</b> 289	579 <b>Uus</b> 293	580 <b>Uuo</b> 293	581 <b>Uuh</b> 289	582 <b>Uus</b> 293
583 <b>Uuh</b> 289	584 <b>Uus</b> 293	585 <b>Uuo</b> 293	586 <b>Uuh</b> 289	587 <b>Uus</b> 293	588 <b>Uuo</b> 293	589 <b>Uuh</b> 289	590 <b>Uus</b> 293	591 <b>Uuo</b> 293	592 <b>Uuh</b>										

## LAB 1A: INTRODUCTION TO CHEMISTRY PRACTICALS

### AIM:

In this introductory practical you would be asked:

1. To identify and name some of the apparatus used in Chemistry Practicals.
2. To measure the volume of a liquid correctly using a measuring cylinder,
3. To measure temperature correctly using a mercury-in-glass thermometer.

### APPARATUS:

- Burette
- Graduated pipette
- Fixed volume pipette
- Retort stand boss and clamp
- Conical flask
- Beaker
- Funnel
- Wash bottle
- Measuring cylinders
- Volumetric flasks
- Glass stirring rod
- Watch glass
- Evaporating basin
- Crucible tongs
- Wash bottle
- Thermometer
- Stopwatch
- Balance
- Test tubes
- Boiling tube
- Test tube rack
- Test tube holder
- Spatula
- Bunsen burner

### MATERIALS:

Distilled water

### PROCEDURE:

- 1) You are provided with various apparatus used in chemistry practicals. The names of these apparatus are listed above. You are to identify, name and draw diagrams/sketches of the apparatus.
- 2) Draw up a list of the names and uses of each apparatus as follows:

Example:

Name of Apparatus	Use
Burette	Measure volumes of liquids

- 3) Take the 100 cm<sup>3</sup> measuring cylinder, place it on the bench top and add water from the wash bottle to give a reading of 45 cm<sup>3</sup> of water measured. You must view the bottom of the water meniscus at eye level.
- 4) Take the thermometer and place it in the water in the measuring cylinder. Record the reading on the thermometer when the mercury level is steady. View the mercury meniscus at eye level.
- 5) Repeat steps 3 and 4 using the 10 cm<sup>3</sup> measuring cylinder (measure 8 cm<sup>3</sup> of water).

#### TREATMENT OF RESULTS:

- 1) Place in order of accuracy, most accurate first, the accuracy of the following apparatus used for measuring liquid volumes: Burette, Beaker, Measuring cylinder and Volumetric flask.
- 2) To what decimal position can you measure accurately, using the thermometer you were provided with?
- 3) Why is accuracy of measurements in the lab so important?
- 4) Why is lab safety so important?
- 5) Distinguish between distilled and deionised water.
- 6) Why do we always view the meniscus at eye level?

## Lab 1B: DETERMINATION OF MELTING POINT

### AIM:

To Determine the Melting-Point of Stearic Acid from its Cooling Curve.

### APPARATUS:

Bunsen burner, tripod, wire gauze, beaker (250 cm<sup>3</sup>), boiling tube, test tube rack, test tube holder, thermometer, stop watch

### MATERIALS:

Stearic acid

### PROCEDURE:

- 1) Set up a water bath by pouring approximately 150 cm<sup>3</sup> of water into a 250cm<sup>3</sup> beaker and then placing it on the tripod and gauze with Bunsen burner under the tripod.
- 2) Place the boiling tube containing the pure powdered stearic acid in the water bath so that the level of water is above the level of the stearic acid in the boiling tube-you may need to top up the water in the beaker.
- 3) Light the Bunsen burner under the water bath and heat until the stearic acid melts.
- 4) When the stearic acid has melted place a thermometer in the boiling tube- see that the thermometer bulb is well covered by the stearic acid without it touching the bottom or sides of the boiling tube.
- 5) Continue heating the liquid stearic acid until it is at about 80 °C.
- 6) Turn off the burner and transfer the boiling tube to the test tube rack
- 7) Take the temperature reading and start the stop clock.
- 8) Take temperature readings every minute until it cools down till you get 3 consecutive readings as the stearic acid cools the thermometer should be GENTLY turned to stir the liquid, but this *must be stopped* as soon as the first crystals of solid appear.
- 9) Note the temperature when solid stearic acid first appears.
- 10) Record your results in a table.
- 11) Plot a graph of Temperature °C (y-axis) against Time in minutes (x-axis).
- 12) At the end of the experiment leave the thermometer in place- do not attempt to remove the thermometer embedded in the solid stearic acid.

### *QUESTIONS:*

- 1) At what temperature did the first crystals of stearic acid appear?
- 2) How can you use your graph to find the freezing point of stearic acid?
- 3) What is the melting point of pure stearic acid?
- 4) If the stearic acid were impure, would the freezing point of stearic acid be higher or lower than the value you obtained in this experiment?
- 5) Explain the shape of your graph.

## Lab 2: Separation of Mixtures

### AIM

To separate mixtures using different separating techniques


### a. Filtration

#### APPARATUS:

Beaker, conical flask, filter paper, funnel

#### MATERIALS:

Sand and water




#### PROCEDURE:

1. Set up your apparatus as seen in the diagram above.
2. Your lecturer will show you how to flute your filter paper.
3. Pour your mixture of sand and water into the funnel slowly. Ensure all the mixture gets into the filter paper.
4. Allow the mixture to pass through the filter paper.
5. Record your observations.

#### QUESTIONS:

1. Define what a mixture is.
2. Which substance remained in the funnel?
3. What is the term given to the substance which remained in the funnel?
4. Why did the substance remain in the funnel?
5. Is the filtrate pure? Explain.
6. Why do you need to flute your filter paper?

## b. Separating funnel



### *APPARATUS:*

Separating funnel and beaker

### *MATERIALS:*

Oil and water


### *PROCEDURE:*

1. Pour the mixture into the separating funnel. Ensure the tap is closed off.
2. Close the top of the separating funnel and shake vigorously.
3. Keep shaking and stop at 20 second intervals and release the air pressure inside the separating funnel by holding the funnel upwards and opening the tap.
4. Close the tap and re-do this about 5 times.
5. Allow to rest until you see a distinct line between the 2 solutions.
6. Gently pour off the lower level solution into a beaker by opening the tap. Stop when the low density solution reaches the tap.
7. For accuracy, discard the area between the lines of separation of the liquids.

### *QUESTIONS:*

1. Give definitions of immiscible and miscible liquids.
2. Give another example of 2 immiscible liquids besides water and oil.
3. Why was the water at the bottom of the mixture?
4. Why did you need to open the separating funnel at different intervals?
5. Why do the solutions not mix?

### c. Chromatography



#### *APPARATUS:*

1000mL Beaker, filter paper, tape or clip, stirring rod or pencil, capillary tube

#### *MATERIALS:*

Dye, salt and water

#### *PROCEDURE:*


1. Cut a piece of filter paper into a 5cm x 10cm rectangle.
2. **Using a pencil**, draw a horizontal line at 1cm from the tip of one end of the filter paper.
3. Place the filter paper into the beaker and allow the bottom of the filter paper to touch the bottom of the beaker and attach the excess outside the top to a stirring rod.
4. Remove the filter paper and pour a 5% salt water solution into the beaker. Fill it to 1cm above the bottom of the beaker.
5. **Using a pencil**, label 2 spots along the horizontal line. Using a dropper, add a drop of each dye onto each spot and wait 5 minutes for it to dry.
6. Place the filter paper into the beaker and ensure that the liquid in the beaker is below the horizontal line. Cover the beaker and leave for 10 minutes or until the solvent reaches about 3 cm from the top of the filter paper.
7. Observe your recordings. Measure the movements and colours seen.

#### *QUESTIONS:*

1. Which is the stationary phase and which is the mobile phase?
2. List the colours which make up each dye.
3. Why did one dye travel further and faster than the other?
4. What is the purpose of the 5% salt water?



#### d. Simple Distillation (Demonstration)



#### *APPARATUS:*

Bunsen burner, simple distillation apparatus, beaker

#### *PROCEDURE:*

The lab technician and lab lecturer will demonstrate the separation of salt and water using the apparatus seen above.

#### *QUESTIONS:*

1. Define the terms “solute” and “solvent”?
2. What is the reason for the thermometer in the apparatus set up?
3. How does the condenser work?
4. What is the purpose of the boiling chips?
5. What is the difference between simple distillation and fractional distillations?
6. Give an example of a commercial use for simple distillation.

## Lab 3: EFFECT OF TEMPERATURE ON THE SOLUBILITY OF A SALT

### AIM:

To investigate how the Solubility of Potassium Nitrate varies with Temperature.

### APPARATUS:

Bunsen burner, tripod, wire gauze, beaker (250 cm<sup>3</sup>), boiling tube, test tube holder, burette and burette stand, thermometer, test tube rack

### MATERIALS:

Pure crystals of Potassium Nitrate (8.0 g in a boiling tube), distilled water

### PROCEDURE:

- 1) Set up a water bath by pouring approximately 150 cm<sup>3</sup> of water into the 250 cm<sup>3</sup> beaker and then placing it on the tripod and gauze with Bunsen burner under the tripod.
- 2) Add 5 cm<sup>3</sup> of distilled water from a burette to the 8.0 g of potassium nitrate in the boiling tube.
- 3) Place the beaker in the water bath so that the level of water is above the level of the potassium nitrate in the boiling tube (you may need to top up the water in the beaker.)
- 4) Light the Bunsen burner under the water bath and heat until all the potassium nitrate dissolves.
- 5) Place the thermometer in the boiling tube, lift the boiling tube out of the water bath and transfer the tube to a test tube rack
- 6) Allow the boiling tube to cool and record the temperature at which crystals *first appear*.
- 7) Add a further 1cm<sup>3</sup> of distilled water from the burette to the boiling tube and repeat steps 3 - 6

**If you overshoot the 1 cm<sup>3</sup> note the precise reading on the burette, you can read the burette to two decimal places, and record the appropriate total volume. E.g. if you added 1.15 cm<sup>3</sup> instead of 1 cm<sup>3</sup> your total volume would be 6.15 cm<sup>3</sup>. Hence your second entry in your *Table of Results*, under the Volume of water/cm<sup>3</sup> column, would be 6.15 and not 6.0**

- 8) Repeat adding further 1cm<sup>3</sup> portions of distilled water to obtain a total of six readings.
- 9) Record results in a table and calculate the solubility in g per 100 g of water at each temperature.
- 10) Plot a graph of solubility, *g per 100 g of water (y-axis)*, against temperature, °C (*x-axis*).
- 11) Draw a smooth curve through the points on the graph.

### RESULTS:

The solubility of a substance is normally expressed in grams of solute per 100 g of water.

You need to calculate the solubility of potassium Nitrate in these units.

Assuming that  $1\text{cm}^3$  of water has a mass of 1 g the calculation is performed as follows:

You started the experiment with 8.0 g of potassium Nitrate and  $5\text{ cm}^3$  of water.

What is the solubility of potassium Nitrate expressed in g per 100 g of water (at the temperature it dissolves)?

**$5\text{ cm}^3$  i.e. 5 g of water dissolves 8 g of potassium Nitrate**

**Hence  $1\text{ cm}^3$  i.e. 1 g of water will dissolve  $8/5$  g of potassium Nitrate**

**Hence  $100\text{ cm}^3$  i.e. 100g of water will dissolve  $8/5 \times 100$  g of potassium Nitrate = 160 g of potassium Nitrate**

i.e.  $(\text{mass of potassium Nitrate} / \text{volume of water}) \times 100 = \text{solubility of potassium Nitrate}$ .

Hence the solubility of potassium Nitrate at the temperature at which crystals first appear = 160 g per 100 g of water

You need to repeat the calculation for each of the other five volumes of water using the actual volume added- see step 8 under PROCEDURE above

### TABLE OF RESULTS:

Volume of water/ $\text{cm}^3$	Mass of potassium Nitrate/ g	Solubility/ g per 100 g of water	Temperature/ $^{\circ}\text{C}$
5.0	8.0	160.0	
	8.0		
	8.0		
	8.0		
	8.0		
	8.0		

### *QUESTIONS:*

- 1) Does the solubility of potassium Nitrate in water increase or decrease as temperature rises?
- 2) Which temperature is chosen as standard when comparing the solubility of different substances?
- 3) Use the graph from the experiment to work out the temperature at which crystals will start to form when a solution of 30 g potassium Nitrate in 100 g water is cooled.
- 4) Use your graph to determine the solubility of potassium Nitrate at 55 °C

## **APPENDIX A**

### **LAB REPORT SHEET FORMAT**

**NAME:**

**PARTNER(S) NAME:**

**DATE:**

**TITLE:**

**AIM:** A brief description of the purpose of the lab.

**APPARATUS:** A list of all the equipment and chemicals that were used in the exercise.

**METHOD:** This is written in past perfect tense, this is a numbered sequence of the steps that were followed to carry out the exercise.

E.g. The addition of iodine to a test tube would be described in the following way:

“The iodine **was added** to the test tube...”

NOT

“They added iodine to the test tube...”

#### **PRECAUTIONS FOR ACCURATE RESULTS AND SAFETY:**

- List of two important precautions which must be followed when using the equipment for the lab exercise.
- For the procedure followed in the virtual lab, suggest two safety measures which must have been observed.

**RESULTS:** Present the results of the experiment (the things you saw or heard). Draw a table wherever possible (title should be written in caps and underlined **above** the table).

**CONCLUSION:** In relation to the aim of the experiment, the deductions made based on the results obtained.

**ANSWERS TO QUESTIONS:** Answer any post lab questions.

## **APPENDIX B**

### **EXAMPLE OF A COMPLETED LAB REPORT SHEET**

**NAME:** Jane Joelly

**PARTNER'S NAME:** Richard Rally

**DATE:** 20.05.2016

**TITLE:** Osmosis

**AIM:** To observe the process of osmosis.

**APPARATUS:** 3 cups, vinegar, corn syrup, water, egg, stop-clock.

**METHOD:**

1. The egg was soaked in vinegar for two days.
2. The egg shell was then removed leaving the soft membrane still holding the yolk and albumen.
3. The egg was then placed in a cup containing enough corn syrup to cover most of it and left for an hour.
4. Any difference in size and appearance was recorded.
5. The egg was then placed in another cup containing just enough fresh water to cover most of it for another two hours.
6. Any difference in size and appearance was recorded.

**PRECAUTIONS FOR ACCURATE RESULTS AND SAFETY:**

Accuracy:

- The stop clock was tested before use to ensure that it was working properly.
- The batteries in the stop clock were changed to ensure that they did not run out during the experiment.

Safety:

- Gloves should have been worn to protect the hands from vinegar.
- Protective glasses were worn to protect the eyes from the vinegar.

**RESULTS:**

After 1 hour in the corn syrup, the egg was smaller than it was at first. There were also wrinkles on the surface.

After 10 minutes in the fresh water, the wrinkles disappeared. The longer the egg stayed in the fresh water, the larger it swelled.

**CONCLUSION:**

In this experiment, osmosis was observed. Osmosis is the movement of water particles from a dilute solution to a more concentrated solution across a semi-permeable membrane. When the egg shell was removed, the membrane left behind was semi-permeable. The albumen was more dilute than the corn syrup, so water moved from the albumen across the semi-permeable membrane and into the syrup. This caused the egg to reduce in size and wrinkle. The albumen was more concentrated than the fresh water however, so water particles moved across the semi-permeable membrane into the egg, making it swell.

**ANSWERS TO QUESTIONS:**

1. Would the same observations be made with an egg with the shell intact?  
No, because the egg shell is not semi-permeable.

2. If the experiment is carried out using the same method but replacing the corn syrup with salt water what observations would be made?

The same observations would be made once the concentration of the salt solution was higher than that of the albumen.

## APPENDIX C

### LAB RUBRIC

This is the Mark Scheme which will be used to score your Lab Reports

SECTION	COMPETENCY	PERFORMANCE RATING			
		Excellent 4	Average 3	Fair 2	Poor 1
<b>Aim</b>	Effective Communication	The aim is clearly stated.	The aim is slightly unclear.	The aim is quite vague.	The aim is not stated.
<b>Apparatus</b>	Knowledge	Apparatus list is complete. Student gives the correct name of all equipment used.	Apparatus list is complete. Student gives the correct name of most equipment used.	Apparatus list is incomplete. Student gives the correct name of most of the equipment they listed.	Apparatus list is incomplete. Student gives the correct name of only a few pieces of the equipment they listed.
<b>Method</b>	Effective Communication	Method is correct, logical, complete and written in appropriate tense.	Method is incorrect, or illogical, or incomplete or written in inappropriate tense.	Method has two - three deficiencies.	Method is incorrect and illogical and incomplete and written in inappropriate tense.
	Technical expertise	Students suggest at least two precautions and two safety measures for each lab which are relevant and accurate.	Students suggest at least one precaution and one safety measure for each lab is relevant and accurate.	Students suggest at least one precaution or one safety measure for each lab which is relevant and accurate.	Students do not suggest any precaution or safety measure which is relevant and accurate.



<b>Results</b>	Effective Communication	Each observation is accurately and completely included.	Most observations are accurately and completely included.	Few observations are accurately and completely included.	No observations are accurately and completely included.
	Technological Expertise	All results are recorded in appropriate format. All calculations, tables and graphs are done correctly.	Most results are recorded in appropriate format or most calculations, tables and graphs are done correctly.	Most results are not recorded in appropriate format or most calculations, tables and graphs are not done correctly.	Most results are not recorded in appropriate format and most calculations, tables and graphs are not done correctly.
<b>Conclusion</b>	Analytical Ability	The chemical principles underlying each observation are fully explained and relevant conclusions are made.	The chemical principles underlying most observations are fully explained and relevant conclusions are made.	The chemical principles underlying few observations are fully explained and relevant conclusions are made.	No observations are fully explained and irrelevant conclusions are made.
<b>Post lab questions/ Post Lab Worksheet</b>	Analytical Ability	All questions answered correctly.	Most questions answered correctly.	Few questions answered correctly.	No questions answered correctly.

## APPENDIX D

### LAB ASSESSMENT SHEET

The marks for each lab you submit will be recorded on one of these sheets by the Lab Lecturer.

<b>i. Student Name:</b>						
<b>ii. Course:</b>						
<b>iii. CRN:</b>						
<b>iv. Lab #:</b>						
<b>v. Date:</b>						
Section	Competency	SCORE				COMMENTS
<b>Aim</b>	Eff. Comm. 1	4	3	2	1	
	<b>TOTAL</b>					
<b>Apparatus</b>	Knowledge	4	3	2	1	
	<b>TOTAL</b>					
<b>Method</b>	Eff. Comm.	4	3	2	1	
	Tech. Exp.	4	3	2	1	
	<b>TOTAL</b>					
<b>Results</b>	Tech. Exp.	4	3	2	1	
	Eff. Comm.	4	3	2	1	
	<b>TOTAL</b>					
<b>Conclusion</b>	Analytical Ability	4	3	2	1	
	<b>TOTAL</b>					
<b>Post Lab questions</b>	Analytical Ability	4	3	2	1	
	<b>TOTAL</b>					
<b>TOTAL LAB MARK /32:</b>						
<b>TOTAL LAB MARK (%):</b>						
<b>LAB LECTURER:</b>						